



JOURNAL OF
APPLIED
CRYSTALLOGRAPHY

Volume 49 (2016)

Supporting information for article:

A visible light excited fluorescence method for imaging protein crystals without added dyes

Tiit Lukk, Richard E. Gillilan, Doletha M. E. Szebenyi and Warren R. Zipfel

Table S1 List of conditions that yielded crystals in the microbatch tray.

Position	Salt	Buffer	Precipitant
A2		0.1M NH ₄ -citrate	40% (w/v) PEG 3000
A3	0.18M tri-ammonium citrate		40% (w/v) PEG 3350
A4/B2	0.02M calcium chloride	0.1M Na-acetate, pH 4.6	30% (v/v) MPD
A5	0.2M magnesium formate		40% (w/v) PEG 3350
A6	0.2M lithium sulfate	0.1M Na-citrate, pH 4.0	40% (w/v) PEG 1000
A7		0.1M Tris-HCl, pH 9.5	40% (w/v) PEG 8000
A10	0.2M sodium formate		40% (w/v) PEG 3350
A12		0.1M NH ₄ -citrate	20% (w/v) PEG 3000
B1	0.18M tri-ammonium citrate		20% (w/v) PEG 3350
B5		0.1M Tris-HCl, pH 9.5	20% (w/v) PEG 8000
B6	0.2M ammonium formate		20% (w/v) PEG 3350
B7	0.2M ammonium chloride		20% (w/v) PEG 3350
B8	0.2M sodium formate		20% (w/v) PEG 3350
B9	0.2M ammonium phosphate	0.1M Tris-HCl, pH 8.5	50% (V/V) MPD
B10	0.2M sodium nitrate		20% (w/v) PEG 3350
B12	0.2M sodium dithionite		20% (w/v) PEG 3350
C1		0.1M Bicine, pH 8.5	20% (w/v) PEG 6000
C2		0.1M HEPES, pH 7.5	10% (w/v) PEG 8000 8% (v/v) ethylene glycol
C3		0.1M Na-cacodylate, pH 6.5	40% (v/v) MPD 5% (w/v) PEG 8000
C4		0.1M Na-citrate, pH 4.0	40% (v/v) ethanol 5% (w/v) PEG 1000
C5		0.1M Na-acetate, pH 4.6	8% (w/v) PEG 4000
C6	0.2M magnesium chloride	0.1M HEPES, pH 7.0	10% (w/v) PEG 8000
C7		0.1M Na-citrate, pH 4.0	20% (w/v) PEG 6000

C10	0.2M tri-sodium citrate		20% (w/v) PEG 3350
C11	0.2M sodium chloride	0.1M Na-citrate, pH 4.0	20% (w/v) PEG 8000
C12	1M lithium chloride	0.1M Na-citrate, pH 4.0	20% (w/v) PEG 6000
D1	0.2M sodium nitrate		20% (w/v) PEG 3350
	0.2M ammonium chloride		
D2		0.1M Hepes, pH 6.5	10% (w/v) PEG 6000
D3	0.8M sodium phosphate	0.1M Hepes, pH 7.5	
	0.8M potassium phosphate		
D5	0.2M Zinc acetate	0.1M Na-acetate, pH 4.6	10% (w/v) PEG 3000
D8		0.1M Bicine, pH 9.0	10% (w/v) PEG 20000
			2% (v/v) 1,4-dioxane
D9	2M ammonium sulfate	0.1M Na-acetate, pH 4.6	
D10			10% (w/v) PEG 1000
			10% (w/v) PEG 8000
E2	0.2M lithium sulfate	0.1M Na-acetate, pH 4.6	30% (w/v) PEG 8000
E3		0.1M Hepes, pH 7.5	70% (v/v) MPD
E4	0.2M magnesium chloride	0.1M Tris-HCl, pH 8.5	20% (w/v) PEG 8000
E5	0.2M lithium sulfate	0.1M Tris-HCl, pH 8.5	40% (v/v) PEG 400
E7	0.17M ammonium sulfate		25.5% (w/v) PEG 4000
			15% (v/v) glycerol
E10	0.04M potassium phosphate		16% (w/v) PEG 8000
			20% (v/v) glycerol
E11	1M tri-sodium citrate	0.1M Na-cacodylate, pH 6.5	
E12	0.2M sodium chloride	0.1M Na-cacodylate, pH 6.5	2M ammonium sulfate
F1	0.2M sodium chloride	0.1M Hepes, pH 7.5	10% (w/v) isopropanol
F3		0.1M Tris-HCl, pH 10.0	40% (v/v) MPD
F6	1M di-ammonium phosphate	0.1M Na-acetate, pH 4.6	
F8		0.1M Bicine, pH 9.0	10% (w/v) PEG 6000
F11	0.05M cesium chloride	0.1M MES, pH 6.5	30% (v/v) Jeffamine M-600

F12	3.2M ammonium sulfate	0.1M Na-citrate, pH 4.0
-----	-----------------------	-------------------------
