



X : parts per Million : 1H

---- PROCESSING PARAMETERS ----
dc_balance(0, FALSE)
sexp(0.2[Hz], 0.0[s])
trapezoid3(0[%], 80[%], 100[%])
zerofill(1)
fft(1, TRUE, TRUE)
machinephase
ppm

Filename = 20201204-PTLC-II-3.jdf
Author = delta
Experiment = single_pulse.ex2
Sample Id = S#561308
Solvent = CHLOROFORM-D
Actual_Start_Time = 24-DEC-2020 20:13:57
Revision_Time = 6-JAN-2021 14:07:44

Comment = single_pulse
Data Format = 1D COMPLEX
Dim Size = 13107
X Domain = 1H
Dim Title = 1H
Dim Units = [ppm]
Dimensions = X
Site = ECX 300
Spectrometer = JNM-ECX300

Field Strength = 7.0586013[T] (300[MHz])
X Acq_Duration = 2.90717696[s]
X Domain = 1H
X_Freq = 300.52965592 [MHz]
X_Offset = 5[ppm]
X_Points = 16384
X_Prescans = 1
X_Resolution = 0.34397631 [Hz]
X_Sweep = 5.63570784 [kHz]
Irr_Domain = 1H
Irr_Freq = 300.52965592 [MHz]
Irr_Offset = 5[ppm]
Tri_Domain = 1H
Tri_Freq = 300.52965592 [MHz]
Tri_Offset = 5[ppm]
Clipped = FALSE
Scans = 8
Total_Scans = 8

Relaxation_Delay = 5[s]
Recvr Gain = 36
Temp Get = 460.0[dC]
X_90_Width = 11.8[us]
X Acq_Time = 2.90717696[s]
X_Angle = 45[deg]
X_Atn = 1.8[dB]
X_Pulse = 5.9[us]
Irr Mode = Off
Tri Mode = Off
Dante Presat = FALSE
Initial Wait = 1[s]
Repetition_Time = 7.90717696[s]